Chemical Indicator

Result Reference Guide

Initial Color

50 mJ/cm²100 mJ/cm²150 mJ/cm²200 mJ/cm²

EN Chemical Indicator

For monitoring UV-C Radiation disinfection processes of surfaces

It has been demonstrated that many microorganisms can be inactivated by UV-C Radiation. DNA and RNA absorb UV-C Radiation, causing changes in their structures and preventing their replication.

The UV-C Radiation dose required to achieve a microbial load reduction depends not only on the particular microbe, but also on additional factors like intensity (power) of the emitting source, distance between emitting source and surface, and time of exposure. Furthermore, surfaces that only receive UV-C Radiation by reflection, get lower doses that may not be sufficient to deactivate the pathogenic microorganisms present. Considering all this, it is of utmost importance to control that all surfaces to be disinfected, receive the minimum UV-C Radiation dose required.

The chemical indicator Chemdye® CD87 has been specially developed to control the efficacy of UV-C Radiation disinfection processes of surfaces. It has been calibrated to experience different color changes depending on the UV-C Radiation dose received, so allowing to quickly and accurately control the efficacy of the disinfection process.

Product Description

The chemical indicator Chemdye® CD87 has been specially developed to control UV-C Radiation disinfection processes of surfaces. It consists of a 15 mm diameter dot of self-adhesive synthetic substrate printed with reactive indicating ink. The reactive indicating ink changes its color depending on the UV-C Radiation dose received.

Instructions for use

1.-Peel off the chemical indicator Chemdye® CD87 from the roll and stick it on the surface which level of disinfection is pretended to be controlled. Use as many indicators as surfaces intended to be verified.

NOTE: it is recommended to distribute chemical indicators Chemdye® CD87 uniformly along the room to be disinfected, placing at least one per square meter. It is recommended to place chemical indicators Chemdye® CD87 at a height not major than 2 m over the floor level.

2.-Run a normal disinfection cycle.

3.-Ended up the process, contrast the final color of each chemical indicator Chemdye® CD87 used against the color chart provided and determine the dose of UV-C Radiation received. The attached documentation chart allows to quickly record and analyze the results obtained.

NOTE: the register of the exposed chemical indicators Chemdye® CD87 must be kept protected from light.

Storage and shelf life

Store in a dry place, protected from light, at a temperature between 10-30 °C, and at a relative humidity between 30-80 %. Do not wet.

Do not store close to disinfectant or sterilizing agents.

The chemical indicator Chemdye® CD87 has an expiration date of 5 years from the date of manufacture if stored under recommended conditions.

Do not use the chemical indicator Chemdye® CD87 after its expiration date.

Disposal

Discard used chemical indicators according to your country's healthcare and safety regulations.

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